

Guidelines for using the NMR facilities at G5 Schrenk Hall

- **Preregister is required** for any access to G5 Schrenk Hall
 - Use google calendar available at the NMR Institute website: <https://chem.mst.edu/research/nmr-institute/>
 - Add your full name as the title of the timeslot you are reserving (don't use abbreviations, email addresses, or titles such as "busy", "unavailable", or "occupied")
 - Once your information is inputted, click under the "edit" tab to add a guest. In the "guest" space, type "400MHz NMR (liquid)" to add to the correct reservation calendar.
- The entrance door to G5 Schrenk Hall must remain closed, and opened only to enter and exit the facilities.
- **Only one person at a time** is allowed in G5 Schrenk Hall (except during Tuesday's weekly maintenance)
- **No chemicals or equipment** (other than NMR samples confined in NMR tubes) may be brought to G5.
- **Do not use liquid or spray disinfecting products on NMR tubes or NMR-tube spinners** (they can leave a residue that might interfere with your measurement or the operation of the spinner, and they may damage the spinners).
- **Use Kimwipes instead.** Discard Kimwipes immediately after use (don't leave them on the lab bench).
- **Clean up and sanitize** all surfaces you have touched or used (including keyboard and mouse covering plastic bags)
- Nitrile gloves must be discarded upon leaving G5 Schrenk Hall.